

# Determination of Cannabinoids in Hemp seed oil using LC/MS/MS with a Novel End-capping C18 Column




Norikazu Nagae<sup>1</sup>, Woo Young Pyo<sup>2</sup>, Etsuko Shearer<sup>3</sup>  
 1. ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka 552-0001 Japan  
 2. Esporalab, 1F 687-4, Gwangsa-dong, Yangju-si, Gyeonggi-do, Korea  
 3. BioNik Inc. 3397-19 Obuchi, Fuji, Shizuoka 417-0801 Japan  
 \*Corresponding author email: nagae@chromanik.co.jp




## Novel End-capping

An End-capping has been recognized to be an important factor for a silica based reversed phase column. In this study, not only bonding with an end-capping reagent but also conversion of silanol groups to siloxane bond by heating were evaluated as an end-capping.




## End-capping with hexamethyldichlorotrisiloxane and TMS on C18 silica



End-capping reagent moves like a *Geometrid caterpillar*, so that a functional group on the tip of the arm can bond with a silanol group which is located anywhere.  
 We named this end-capping method as Sunniest double end-capping.

## Stability under acidic pH condition




Column: SunShell C18 2.6 μm, 50 x 2.1 mm  
 Mobile phase: A) 0.1% trifluoroacetic acid pH 2.0  
 B) Acetonitrile  
 Gradient program: 

Time (min)	0	3	3.1	5
% B	10	90	10	10

  
 Flow rate: 0.5 mL/min  
 Temperature: 40 °C  
 Detection: UV@270nm  
 Injection volume: 0.5 μL  
 Sample: 1 = Benzylamin (0.5 mg/mL), 2 = Ketoprofen (0.04 mg/mL), 3 = Indomethacin (0.05 mg/mL)

## Bleeding test using LC/MS




Column size: 50 x 2.1 mm  
 Mobile phase: A) 0.1% acetic acid  
 B) Acetonitrile  
 Gradient program: 




Time (min)	0	1	5	7
% B	5	5	10	10

  
 Flow rate: 0.4 mL/min  
 Temperature: 40 °C  
 MS: ABI API-4000  
 Ionization: Turboionspray (cation)  
 Measurement mode: Q1 Scan m/z 100-1000

## Separation of 9 kinds of cannabinoids



## Measurement of CBD and THC in hemp seed oil



Column: Thunder C18 2.7 μm, 150 x 2.1 mm  
 Mobile phase: A) 0.1% Formic acid in water  
 B) 0.1% Formic acid in acetonitrile

Time (min)	0	0.5	2.5	6	6.1	10
%B	30	30	80	80	30	30

Flow rate: 0.4 mL/min  
 Temperature: 40 °C  
 Injection Volume: 5 μL  
 Peaks: A = Cannabidiol (CBD, 1000 μg/mL)  
 B = δ-9-Tetrahydrocannabinol (THC, 10 μg/mL)  
 Sample: Hemp seed oil

Table 2. CBD / THC MRM Transitions

Compound Name	R.T.	Precursor ion	Product ion	CE (V)	Polarity
CBD	3.49	315	193	16.0	Positive
CBD	3.49	315	259	16.9	Positive
THC	3.92	315	193	22.6	Positive
THC	3.92	315	259	22.5	Positive

Table 3. Measurement result of hemp seed oil

	CBD	THC
Average (ppb)	44.2	3.0
Dilution	20	20
Calculation of concentration (ppb)	884 (n=10)	60

**Conclusion**

- Hexamethyldichlorotrisiloxane and trimethylchlorosilane were used as end-capping reagents. Silanol groups were changed to siloxane bonding at an elevated temperature. As a result, stability increased drastically.
- Nine kinds of cannabinoids were successfully separated using the C18 with novel end-capping.
- The concentrations of CBD and THC in hemp seed oil determined using LC/MS/MS were 884 ppb and 60 ppb, respectively.