A Novel Bonding Technique Using a Polyfunctional Silyl-Reagent for Reversed-Phase Liquid Chromatography--- A NEW APPROCH !!!! II

Norikazu Nagae¹⁾, Tomoyasu Tsukamoto¹⁾, Anders Grahn²⁾, Ralf Jutvik²⁾ and Vinay D. Gaitonde²⁾

1) ChromaNik Technologies Inc., 6-3-1 Namiyoke, Minato-ku, Osaka, Japan

Novel C18 silyl-reagent (HMODTS)

Evaluation of end-capping

Comparison of amitriptyline peak I

N(4)=11,000

Patent pending

Step 2

na nata nata na

Column size: 150 X 4.6 mm

2) Biotech AB, Box 133, 439 23, Onsala, Sweden

Step 1

CH.OH. pH7.5. 40 °C

ChromaNik Technologies Inc

Final TMS

:: 5 μm ter: 12n a: 340 m

Bonding state of HMODTS on silica

Evaluation of end-capping

Comparison of amitriptyline peak II

Column size: 150 X 4.6 mm

on the

CH.CN. pH7.0.40 °C

N(4)=14,00 TF(4)=1.16

Abstract

Reversed- phase LC columns and Sorbent development has experienced chromatographic performance improvement in terms of ULTRA pure silica, New End capping Reagents and chemistries, along with improvement in Bonding technology, innovative hybrid silica particle and high speed fused core Sorbents as well as sub 2 micron high efficiency sorbents and columns. These new innovative contributions are well recognized and accepted in HPLC work practices.

In regards to Bonding chemistries, most of reversed phase silica materials are monomerically or polymerically bonded with alkyl chain, then end-capped with trimethylsilane or hexamehyltrisiloxane etc. In this study, polyfunctional silyl-reagent was synthesized with

octadecyltrimethoxylsilane and hexamethyldichlorotrisiloxan. This reagent is called hexamethyloctadecyltetrasiloxane (HMODTS - C18). We are able to bond the Chromatographic Silica particles with this reagent and finally end-cap with trimethylchlorosilane .

The resultant Sorbent was evaluated for chromatographic performance for resolving acidic and basic compounds in various Mobile phase compositions and Buffers. Stability of this phase was evaluated under both acidic and basic pH at high temperature. This phase showed symmetrical peaks of both acidic and basic compounds such as Formic acid and Amitriptyline. Interestingly enough we observed a symmetric peak of Amitriptyline in Mobile phase consisting of Acetonitrile and Ammonium acetate compositions. Majority of the C18 columns showed a long tailing peak of Amitriptyline under similar conditions. We observed that the Column life was more than 500 hours from pH 1.5 to pH 10 at 50 degree Celsius. A novel bonding technique using a polyfunctional silylreagent(HMODTS - C18) could make effect of residual silanol groups the least and hence offer a tailing free chromatographic performance as well as endurance at elevated temperatures.

