


# Some Hint on How to Make a Standard UHPLC Column with + 300 000 Theoretical Plates/Meter



WWW.BIOTECH-AB.COM

Norikazu Nagae<sup>1</sup>, Tomoyasu Tsukamoto<sup>1</sup>, Henrik Svennberg<sup>2</sup>  
<sup>1</sup> ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan  
<sup>2</sup> Biotech AB, Innovation, Onsala, Sweden


**ChromaNik**  
 ChromaNik Technologies Inc.  
 www.chromanik.co.jp



The race to achieve the highest number of theoretical plate counts has been going for many years already, even though it could and should be argued that a good separation is mainly governed by good selectivity, a reasonable retention and then high column efficiency. However, the race has still been going, and this presentation will show how + 300 000 theoretical plates/meter can be achieved in a standard column for UHPLC. The importance of a good bonding and end capping technique will be discussed and hints about how this is made is revealed. The importance of base particle choice with a good heat transfer capacity will be shown and the importance of small extra column peak broadening will be emphasized. The + 300 000 theoretical plates/meter is shown in two different columns lengths with the inner diameter of 2.1 mm, 150 and 50 mm long respectively.


## Effect of end-capping

### End-capping with hexamethyltrisiloxane and TMS on C18 silica



End-capping reagent moves like a *Geometrid caterpillar*, so that a functional group on the tip of the arm can bond with a silanol group which is located anywhere. We named this end-capping method **Sunniest end-capping**.


### Another end-capping with heating on C18 silica, reduce of silanol groups



A basic compound shows no tailing on Sunrise C18-SAC because hydrated silanol groups don't make a basic compound tailing as well as silica column on HILIC mode shows no tailing for a basic compound.


- Non hydrated silanol group by influence of hydrophobicity of alkyl groups
- Hydrated silanol group without influence of alkyl groups

### Comparison of 2 kinds of end-capping



### Carbon loading and retention factor (reaction at two different temperatures)

End-capping reaction temperature	X °C (more than 200 °C)	(X-40) °C
Carbon loading of only C18	7.0%	7.0%
Carbon loading after end-capping	7.3% (not cut off C18 chain by heat)	7.7% (not cut off C18 chain by heat)
Silanol activity control	Yes	No
Retention factor (k)	10.6	9.5




### Relationship between retention factor and carbon loading

Column	Retention factor <sup>a</sup>	Carbon loading (%)	Specific surface area <sup>b</sup> (m <sup>2</sup> /g)
SunShell C18, 2.6 μm	10.4	7.3	125
Ascentis Express C18, 2.7 μm	9.7	8.0	133
PoroShell C18 EC, 2.7 μm	9.0	8.5	135
Cortex C18, 2.7 μm	7.7	7.3	113
Accucore C18, 2.6 μm	7.4	8.8	130
Kinetex C18, 2.6 μm	5.4	4.9	102

<sup>a</sup> Retention factor of acetonitrile, mobile phase, acetonitrile/10mM ammonium acetate pH6.8 = 40/60; flow rate, 1.0 mL/min; temp.: 40 °C.  
<sup>b</sup> Measured after C18 packing material was pulled out of a column.  
<sup>c</sup> Measured using C18 materials sintered at 600 degree Celsius for 8 hours.

Both Sunniest end-capping and Silanol Activity Control

### Comparison of amitriptyline peak




Column dimension, 150 x 4.6 mm; mobile phase, acetonitrile/10mM ammonium acetate pH6.8 = 40/60; flow rate, 1.0 mL/min; temp.: 40 °C; sample, 1-uracil, 2-propranolol, 3-nortriptyline, 4-amitriptyline.

SunShell C18 showed the highest theoretical plate (N) and the lowest tailing factor. Highly-effective end-capping including heat treatment (silanol activity control) reduced tailing for a basic compound. As a consequence, high efficiency for a reversed column was achieved by highly-effective end-capping.

## Effect of frictional heating


### Frictional heating phenomenon



Column, core shell C18 2.6 μm, 150 x 4.6 mm; mobile phase, methanol, flow rate, 5 mL/min; column back pressure, 70 MPa; temperature, 23 °C (room temperature).


Temperature rose 14 °C at the outlet of the column by effect of frictional heating.

### Efficiency loss by thermal friction\*




Column dimension, 100 x 2.1 mm; mobile phase, CH<sub>3</sub>CN/H<sub>2</sub>O=60/40; temperature, 40 °C; sample, acenaphthene; UHPLC instrument, X-LC with 0.1 mm i.d. connecting tube (JASCO).


### Efficiency and back pressure



### Peak shape and efficiency of pyridine and amitriptyline



Column dimension, 100 x 2.1 mm; mobile phase, methanol/water = 30/70; flow rate, 0.2 mL/min; temperature, 40 °C; sample, 1-uracil, 2-pyridine, 3-phenol; UHPLC instrument, X-LC with 0.1 mm i.d. connecting tube (JASCO).



Column dimension, 100 x 2.1 mm; mobile phase, acetonitrile/10mM ammonium acetate pH6.8 = 40/60; flow rate, 0.3 mL/min; temperature, 40 °C; sample, 1-uracil, 2-pyridine, 3-phenol; UHPLC instrument, X-LC with 0.1 mm i.d. connecting tube (JASCO).

Influence of frictional heating is most difficult to undergo for a core shell silica. When using a fine particle such as 2 μm or smaller than 2 μm C18 packings under high pressure condition, a core shell C18 with dense end-capping showed the highest efficiency for not only neutral compounds but also basic compounds.

## Effect of extra-column volume


### Comparison of an inner diameter of tubing in the UHPLC instrument

Connecting tube used	Injector to column	Column to flow cell	Tube volume
X-LC with SUS	SUS 0.1 mm i.d. x 300 mm	PeekSil 0.1 mm i.d. x 200 mm	3.93 μL
X-LC with MarvelX	MarvelX 0.075 mm i.d. x 350 mm	MarvelX 0.075 mm i.d. x 150 mm	2.21 μL

\*X-LC, is a UHPLC instrument manufactured by Jasco. \*MarvelX is a connecting tube manufactured by GEEK.

### SunShell C18 2.6 μm, 50 x 2.1 mm

#### X-LC with SUS




Measurement condition  
 Column, SunShell C18, 2.6 μm 50 x 2.1 mm; mobile phase, CH<sub>3</sub>CN/H<sub>2</sub>O=60/40; flow rate, 0.6 mL/min; temperature, RT; detection, UV@250 nm; injection volume, 0.4 μL; sample, 1-uracil, 2-ethylbenzene, 3-acenaphthene, 4-butylbenzene.

Peak No.	SUS	Marvel X
Efficiency	1 1208 3593 197% up	2 7720 12625 64% up
Tailing factor	3 1.9589 1.5153 12% up	4 1.3936 1.4733 6% up
Peak width, W <sub>0.5</sub> (min)	1 0.0094 0.0051	2 0.0107 0.0083
Peak width, W <sub>0.5</sub> (min)	3 0.0201 0.0194	4 0.0297 0.0295

### SunShell C18 2.6 μm, 150 x 4.6 mm

[[An extremely low dead volume column hardware manufactured by Tomoe Works was used.]]

#### X-LC with SUS




Measurement condition  
 Column, SunShell C18, 2.6 μm 150 x 4.6 mm; mobile phase, CH<sub>3</sub>CN/H<sub>2</sub>O=70/30; flow rate, 1.8 mL/min; temperature, RT; detection, UV@250 nm; injection volume, 0.4 μL; sample, 1-uracil, 2-ethylbenzene, 3-acenaphthene, 4-butylbenzene.

Peak No.	SUS	Marvel X
Efficiency	1 41900 48932 17% up	2 42934 44037 2.7% up
Tailing factor	3 1.3889 1.4089 4.9% up	4 1.3825 1.3947 4.8% up
Peak width, W <sub>0.5</sub> (min)	1 0.0074 0.0070	2 0.0172 0.0170
Peak width, W <sub>0.5</sub> (min)	3 0.0395 0.0390	4 0.0542 0.0533

### SunShell C18 2.6 μm, 150 x 2.1 mm

#### X-LC with SUS



Measurement condition  
 Column, SunShell C18, 2.6 μm 150 x 2.1 mm; mobile phase, CH<sub>3</sub>CN/H<sub>2</sub>O=70/30; flow rate, 0.4 mL/min; temperature, 25 °C; pressure, 64 MPa; detection, UV@250 nm; injection volume, 0.4 μL; sample, 1-uracil, 2-ethylbenzene, 3-acenaphthene, 4-butylbenzene.

Peak No.	SUS	Marvel X
Efficiency	1 41900 48932 17% up	2 42934 44037 2.7% up
Tailing factor	3 1.3889 1.4089 4.9% up	4 1.3825 1.3947 4.8% up
Peak width, W <sub>0.5</sub> (min)	1 0.0074 0.0070	2 0.0172 0.0170
Peak width, W <sub>0.5</sub> (min)	3 0.0395 0.0390	4 0.0542 0.0533

## Conclusion

- End-capping done by both bonding of a silylation agent and conversion from silanol groups to siloxane bond led to not only high efficiency for basic compounds but also increase in retention.
- Serious losses in column efficiency was caused by frictional heating effects when very fine particles were packed into a column and it was operated at high back pressure. A core shell particle showed the least losses in column efficiency, while a totally porous hybrid particle showed large losses.
- Effect of extra-column volume is important for a column packed with very fine particles. When 0.075 mm i.d. connecting tubes were used, SunShell C18 2.6 μm column showed + 300,000 theoretical plates/meter.
- In case of a low dead volume column sized 4.6 mm i.d. and 150 mm length, uracil peak showed + 300,000 theoretical plates/meter even if a particle size was 2.6 μm.
- In case of 150 mm length column, efficiency was + 300,000 theoretical plates/meter even if 0.1 mm i.d. connecting tubes were used.