

## カンナビノイドCBDとTHCの分離

### Cannabidiol and $\delta$ -9-Tetrahydrocannabinol

SunShell C18 2.6  $\mu$ m, 150 x 2.1 mm i.d.




Fig 1. (A), (B) : 1, 2.5, 5, 10, 50, 100 ppb STD Overlay





Fig 2. (A), (B) : TIC (Total Ion Chromatogram) / Spectra

#### CBD

$y = 8210x + 1120$   
 $R^2 = 0.99954$   
 RSD RF 14.48  
 Creation Date 2021-06-01 09:50:26  
 Internal Standard  
 Regression Linear  
 Weighting None  
 Origin IGNOREZERO  
 Signal AREA  
 Min Conc 1.0000 ppb  
 Max Conc 100.00 ppb



#### THC

$y = 7179x + 4902$   
 $R^2 = 0.99984$   
 RSD RF 7.989  
 Creation Date 2021-06-01 12:54:25  
 Internal Standard  
 Regression Linear  
 Weighting None  
 Origin IGNOREZERO  
 Signal AREA  
 Min Conc 1.0000 ppb  
 Max Conc 100.00 ppb






Fig 3. Verification of Calibration Curve (1 ~ 100ppb)



Cannabidiol (CBD)



$\delta$ -9-Tetrahydrocannabinol (THC)

#### CBD / THC MRM Transitions

Compound Name	R.T	Precursor ion	Product ion	CE (V)	Polarity
CBD	3.49	315	193	16.0	Positive
CBD	3.49	315	259	16.9	Positive
THC	3.92	315	193	22.6	Positive
THC	3.92	315	259	22.5	Positive

#### Measurement result of hempseed oil

	CBD	THC
Average (ppb)	44.2	3.0
Dilution	20	20
Calculation of concentration (ppb)	884 (n=10)	60

Column: SunShell C18 2.6  $\mu$ m, 150 x 2.1 mm

Mobile phase:

A) 0.1% Formic acid in water

B) 0.1% Formic acid in acetonitrile

Time (min)	0	0.5	2.5	6	6.1	10
%B	30	30	80	80	30	30

Flow rate: 0.4 mL / min

Temperature: 40  $^{\circ}$ C

Injection Volume: 5  $\mu$ L

Peaks: A = Cannabidiol (CBD, 1000  $\mu$ g/mL)

B =  $\delta$ -9-Tetrahydrocannabinol (THC, 10  $\mu$ g/mL),

Sample, Hempseed oil

MS (ChroZen LC-TQD) conditions

Source type: HESI

Spray voltage: 4500 V (Positive)

Cone temperature: 250  $^{\circ}$ C

Cone gas pressure: 20 psi

Heated probe temperature: 300  $^{\circ}$ C

Probe gas pressure: 40 psi

Nebulizer gas pressure: 40 psi

Active exhaust: ON

Courtesy of Young In and Esporalab