

# Evaluation of Durability of Hybrid Silica with Ethylene Chains



**BioNik Inc.**

www.bionikinc.com

N. Nagae<sup>a</sup>, T. Tsukamoto<sup>a</sup>, R. Koyama<sup>a</sup>, E. Shearer<sup>b</sup>

[a] ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan

[b] BioNik Inc. Obuchi, Fuji, Shizuoka Japan



www.chromanik.co.jp



Silica-based reversed-phase columns have been widely used since the 1970s. Silica has the advantage of being mechanically strong and able to withstand high pressures, but conventional silica C18 has limited use under alkaline conditions and is only durable up to pH 8. In the 2000s, hybrid silica with ethylene chains bonded to the inside or surface of silica was developed, and patent applications were filed for these technologies. Currently, these technologies are freely available, and it is expected that these hybrid silica packings will be available on the market.

In this study, we evaluated the durability of so-called hybrid silica with ethylene chains bonded to the inside and surface of silica, and hybrid C18 bonded with C18 and end-capped. The amount of hybrid silica dissolved was measured when an alkaline mobile phase of pH 11.5 was passed through the column, and the hybrid C18 bonded with C18 and end-capped was also evaluated in the same way. The durability of hybrid silica under alkaline conditions greatly differed depending on the method of introducing and bonding ethylene chains to silica, with a difference of 100 to 1000 times the amount of dissolved silica. The difference in stability of hybrid C18 under alkaline conditions was smaller than that of hybrid silica alone, suggesting that the effect of end-capping on durability is greater than that of hybrid silica.

## 2 types of hybrid silica gel



### 1. Hybrid ①

Ethylene cross-linked silica gel



### 2. Hybrid ②


Ethylene chains bonded to the silica gel surface



## Hybrid silica gels used in this study

Name	Specific surface area	Pore diameter	Carbon loading	reagent
Hybrid ① (6.4%)	180 m²/g	15 nm	6.4%	
Hybrid ② (2.3%)	340 m²/g	12 nm	2.3%	Silica gel +
Hybrid ② (1.2%)	340 m²/g	12 nm	1.2%	Silica gel +
Silica with end-capping (4.4%)	340 m²/g	12 nm	4.4%	Silica gel +

## Comparison of durability of hybrid silica gel under alkaline conditions



### Durability test conditions (pH 11.5)

Column dimension: 100 x 2.1 mm  
Mobile phase: 50 mM potassium phosphate pH 11.5

Flow rate: 0.1 mL/min

Temperature: 40 °C

Elution time: 17 h

(The carbon loading of the ethylene chain is shown in parentheses.)

Hybrid ① (6.4%) } Full hybrid


Company A (7%) } Full or partial hybrid

Company B (7.3%) } Full or partial hybrid

Hybrid ② (2.3%) } Partial hybrid

Hybrid ② (1.2%) } Partial hybrid

Silica with only end-capping (4.4%) } No hybrid



### Durability test conditions (pH 11.5)

Column dimension: 100 x 2.1 mm

Mobile phase: 50 mM potassium phosphate pH 11.5

Flow rate: 0.1 mL/min

Temperature: 40 °C and 60 °C



Elution time: 17 h

- When comparing the amount (speed) of silica dissolution when an alkaline mobile phase was passed through hybrid silica, Hybrid ② showed a value 100 to 1000 times higher than Hybrid ①.
- The amount of silica C18 dissolved in an alkaline mobile phase increases exponentially over time, and this was also the case for hybrid silica (Hybrid ①).
- When the temperature difference was 20°C, the amount of silica dissolved was about 10 times different, and a 10°C temperature increase caused the amount of silica dissolved to increase by about three times.

## C18 bonding and end-capping

## Durability under acidic and basic pH conditions

## Conclusions



- ✓ It was confirmed that the durability of hybrid silica under alkaline conditions varies greatly depending on the method of bonding the ethylene chains.
- ✓ It was suggested that the durability of the C18 packing after C18 bonding and end-capping is greatly affected not only by the durability of the hybrid silica, but also by the surface treatment, especially end-capping.
- ✓ C18 using hybrid silica with ethylene chains bonded within the silica skeleton was found to have sufficient durability under mobile phase conditions from pH 1 to pH 12.