

Evaluation of Hybrid Silica C18 End-capped with Bidentate Silylating Reagent for HPLC



BioNik Inc.

www.bionikinc.com


N. Nagae^a, T. Tsukamoto^a, R. Koyama^a, E. Shearer^b

[a] ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan


[b] BioNik Inc. Obuchi, Fuji, Shizuoka Japan

ChromaNik
ChromaNik Technologies Inc.
www.chromanik.co.jp

Comparison between Company W hybrid C18 and SunBridge C18 with Ethylene cross-linked silica gel And 1,2-Bis(chlorodimethylsilyl)ethane end-capping



Evaluation of hydrogen bonding, hydrophobicity and steric selectivity



As a comparison, Hybrid C18 manufactured by W Company, which is ethylene crosslinked silica gel C18 in which an ethylene chain is incorporated into the silica skeleton, was used.

Column: Company W Hybrid C18 3.5 µm, 5 µm,

SunBridge C18 5 µm,


Column demensions: 150 x 4.6 mm

Mobile phase: CH₃OH/H₂O=75/25

Flow rate: 1.0 mL/min


Temperature: 40 °C

Sample: 1 = Uracil, 2 = Caffeine, 3 = Phenol, 4 = Butylbenzene,



	Hydrogen bonding α(Caffeine/Phenol)	Hydrophobicity α(Amylbenzene/ Butylbenzene)	Steric selectivity α(Triphenylene/ o-Terphenyl)	Specific surface area (m ² /g)	Carbon loading (%)
Company W Hybrid C18 3.5 µm	0.38	1.54	1.33	185	17.4%
Company W Hybrid C18 5 µm	0.39	1.52	1.33	188	18.0%
SunBridge C18 5 µm	0.43	1.53	1.31	190	16.0%


Comparison of peaks for pyridine, oxine, and nortriptyline



Company W Hybrid C18 showed severe tailing of metal-chelating, oxine.

Compared to Company W Hybrid C18, SunBridge C18 showed a higher number of plates and less tailing of peaks for metal-chelating and basic compounds.

Comparison of Stability under pH 11.5 and pH 1.0 conditions



Durable test condition (pH 11.5)
Column dimension: 150 x 4.6 mm
Mobile phase: Methanol/50 mM potassium phosphate pH 11.5=10/90
Flow rate: 1.0 mL/min
Temperature: 40 °C


Measurement condition
Mobile phase: Acetonitrile/water=70/30
Flow rate: 1.0 mL/min
Temperature: 40 °C
Sample: 1 = Uracil, 2 = Butylbenzene

Durable test condition (pH 1.0)
Column dimension: 50 x 2.1 mm
Mobile phase: Acetonitrile/1.0% TFA pH1=10/90
Flow rate: 0.4 mL/min
Temperature: 80 °C

Measurement condition
Same as above

SunBridge C18 showed almost the same or higher stability under pH 11.5 and pH 1.0 conditions as Company W Hybrid C18, which shows an order of magnitude higher stability than other companies' hybrid C18.

Reproducibility in retention under 100% aqueous conditions



Conclusions

- ✓ SunBridge C18 was developed using Ethylene cross-linked silica gel and 1,2-Bis(chlorodimethylsilyl)ethane end-capping reagent.
- ✓ As a comparison, Hybrid C18 manufactured by W Company, which is ethylene cross-linked silica gel C18 in which an ethylene chain is incorporated into the silica skeleton, was used.
- ✓ Company W Hybrid C18 showed severe tailing of metal-chelating, oxine. Compared to Company W Hybrid C18, SunBridge C18 showed a higher number of plates and less tailing of peaks for metal-chelating and basic compounds.
- ✓ SunBridge C18 showed almost the same or higher stability under pH 11.5 and pH 1.0 conditions as Company W Hybrid C18, which shows an order of magnitude higher stability than other companies' hybrid C18.
- ✓ SunBridge C18 showed an excellent reproducibility in retention factor under 100% aqueous conditions while Company W Hybrid C18 showed a poor reproducibility.