


Analysis of PFAS in Tap Water Using a Pentafluorophenyl Column



Tadashi Kitta¹, Hirotake Takahashi¹,
Norikazu Nagae², Ryuji Koyama² and Tomoyasu Tsukamoto²
¹Japan Food Inspection Corporation, Heiwajima, Ohta-ku, Tokyo, Japan
²ChromaNik Technologies Inc., Namiyoke, Minato-ku, Osaka, Japan




Perfluoroalkyl and polyfluoroalkyl substances (PFAS) are known to pose health risks due to their bioaccumulative nature and environmental persistence. In April 2026, legal restrictions setting a combined limit of 50 ng/L for PFOS and PFOA in tap and mineral water will take effect in Japan. Furthermore, the European Union has already established maximum levels for four PFAS congeners in seafood, and similar regulations are anticipated in the United States.

In this study, we developed a pentafluorophenyl (PFP) column consisting of an ethylene cross-linked hybrid silica gel and a PFP stationary phase. This column is capable of simultaneously analyzing 19 PFAS compounds, ranging from short-chain to long-chain substances. Its applicability to drinking water analysis was evaluated and compared with that of a conventional C18 column. The PFP column exhibited not only greater sensitivity for PFBA but also superior separation of PFHxS isomers compared to the C18 column.

PFAS


Per- and polyfluoroalkyl substances (PFAS) are defined as "fluorinated substances that contain at least one fully fluorinated methyl or methylene carbon (with no H/Cl/Br/I atoms attached). As of January 2025, there are 14,735 PFAS with clear chemical structures and 1,915 PFAS with unclear structural information such as polymers and complex mixtures. PFAS are heat-resistant, water-repellent, and oil-repellent, and are used in many industrial products and daily necessities. On the other hand, they are called "Forever Chemicals" because they are difficult to decompose in the environment and higher bioaccumulative. The main routes of human exposure are oral ingestion of drinking water and food, inhalation of product dust, and skin absorption through contact with products (such as cosmetics).



Hybrid silica PFP and C18


1. Hybrid silica

Ethylene cross-linked silica gel




2. Stationary phase

PFP (SunBridge PFP-R)



C18 (SunBridge C18)



+ End-capping

Analysis issues

- When a conventional C18 column is used,
- 1) lack sensitivity for short-chain PFAS.
 - 2) insufficient separation of long-chain PFAS isomers.
- A method that can comprehensively and sensitively analyze a variety of PFASs is needed.


Objectives

Aim to establish a comprehensive and highly sensitive method for analyzing trace PFAS in tap water using PFP columns.

Measurement conditions

HPLC	Ultimate 3000(Thermo Fisher Scientific)
Analytical column	SunBridge PFP-R 3 μm, SunBridge C18 3 μm, 100 x 2.1 mm ID
Delay column	SunBridge Delay C18 3 μm, 50 x 3.0 mm ID
Mobile phase	A) Methanol B) 2 mM Ammonium acetate
Gradient	Time (min) %B 0.0 30 1.0 30 21.0 97 24.0 97 24.1 30 28.0 30
Flow rate	0.3 mL/min
Column temperature	40 °C
Injection volume	4 μL
MS	Q Exactive Focus (Thermo Fisher Scientific)
Ion mode	ESI (-)
Resolution	70000 (at m/z 200)


Sample pretreatment



Analytes

PFOS, PFOA, PFBS, PFHxS, PFBA, PFPeA, PFHxA, PFHpA, PFNA, HFPO-DA

Comparison of isomer separation




Comparison of PFBA peak intensity between PFP column and C18 column


Column	Peak intensity		
	100% MeOH	30% MeOH +70% ammonium formate	100% MeCN
SunBridge PFP-R	442,000	173,000	336,000
SunBridge C18	277,000	110,000	218,000

Isocratic separation using 100% MeOH, 100% MeCN and 30% MeOH + 70% ammonium formate. Sample: PFBA (5 ng/mL), injection volume: 4 μL

Separation of 19 kinds of PFASs in standards and water sample



Comparison of separation of 10 kinds of PFASs



Conclusions

- ✓ SunBridge PFP-R exhibited not only greater sensitivity for PFBA but also superior separation of PFHxS isomers compared to SunBridge C18.
- ✓ SunBridge PFP-R was able to completely separate 19 kinds of PFASs although the C18 was unable to adequately separate PFNA and PFOS.
- ✓ SunBridge PFP-R is made of hybrid silica gel, so it is expected to have a longer column life than conventional PFP columns.