

SunBridge



**PFAS -
analysis
compatible**


PFP-R



SunBridge PFP-R (3 μm)

New

**SunBridge Innovation
for
High-Stability PFP!**



Ideal for short- and long-chain PFAS analysis

100%MeOH

PFP-R


C18

PFBA

PFPeA

PFOA

PFOS




- Column focusing — peak shape maintained
- Temperature screening — impurity separation
- LOD $\leq 0.05 \mu\text{g/kg}$ for all target PFAS
- Clear separation — Branched and linear PFAS

For more information on PFAS analysis:

https://chromanik.co.jp/info/wp-content/uploads/2025/07/hplc2025_poster_d.pdf

The Redesigned **PFP** for **R**eversed-phase



Ultra Hybrid Technology




PFP-R: Redesigned **PFP** for **R**eversed phase chromatography

Features of the PFP-R Stationary Phase


- Dipole–dipole interactions based on the strong electron-withdrawing property of fluorine
- π/π donor–acceptor interactions (strong recognition of aromatic rings)
- Tuned to provide moderate (rather than strong) retention for cations
- PFP designed with emphasis on stability as a reversed-phase stationary phase
- Note: The so-called “HILIC-like behavior” often attributed to PFP columns (see figure below) does not occur.



Only for basic compounds; driven by silanol, not the PFP phase¹⁾.
Traditional PFP shows a U-shape, but PFP-R does not.



Cation retention test



Column :

SunBridge PFP-R 3 μ m Batch#1

SunBridge PFP-R 3 μ m Batch#2

Sunniest PFP 5 μ m, SunBridge C18 3 μ m

Column dimension: 150 x 4.6 mm

Mobile phase: (A) 0.2% formic acid (pH 2.6)

(B) 100 mM Ammonium formate buffer (pH 3.6)

(C) 100 mM 100 mM Ammonium formate (pH 6.3)

Flow rate: 1.0 ml/min, Temperature: 25 $^{\circ}$ C,


Detection: UV@260 nm

Sample: (a)=L-Adrenaline, (d)=L-DOPA


Cation Retention indicator: $\alpha(a/d)$

- Cation retention was lower than that of traditional PFP, with minimal batch-to-batch variation.

Chelation test




Formic acid test



- Formic acid and Oxine, despite strong adsorption and peak tailing tendencies, showed good peak shapes.

Alkaline durability test



Durability test condition

Column: SunBridge PFP-R 3 μ m, 50 x 2.1 mm

Traditional PFP(PFPP) 5 μ m, 50 x 2.0 mm

Mobile phase: 20 mM H₃PO₄ buffer(pH 8.0)

Flow rate: 0.2 mL/min

Temperature: 40 $^{\circ}$ C

Sample: butylbenzene

Theoretical plate measurement condition

Column: SunBridge PFP-R 3 μ m, 50 x 2.1 mm

Traditional PFP(PFPP) 5 μ m, 50 x 2.0 mm

Mobile phase: CH₃OH/H₂O = 50/50

Flow rate: 0.2 mL/min

Temperature: 40 $^{\circ}$ C

Detection: UV@250nm

Sample: butylbenzene

- ✕ PFPP exceeded pressure limits, making it impossible to continue the test within 48 hours.
- ⊙ SunBridge PFP-R exhibited minimal loss of theoretical plates after 300 hours of flow.*

* Accelerated test confirmed substrate stability: no loss of plates, but reduced retention time.
 As with other PFP columns, use under acidic conditions is recommended to suppress silanol activity.

Organofluorine selectivity


SunBridge PFP-R (Batch#1)

$\alpha(\text{FB})$

=1.10

$\alpha(\text{TFT})$

=1.34





Sunniest PFP

$\alpha(\text{FB})$


=1.14

$\alpha(\text{TFT})$

=1.49



PFP-R



Column :

SunBridge PFP-R 3 μm Batch#1

SunBridge PFP-R 3 μm Batch#2

Sunniest PFP 5 μm

SunBridge C18 3 μm

Column dimension: 150 x 4.6 mm

Mobile phase: Methanol/Water(60 : 40)

Flow rate: 1.0 ml/min, Temperature: 40 $^{\circ}\text{C}$,

Detection: UV@250 nm

Sample: 1=Uracil(t0), 2=Benzene,

3=Fluorobenzene (FB),

4=Toluene, 5=Bromobenzene,

6= α,α,α -Trifluorotoluene (TFT)

Organofluorine Selectivity Indicator-1 : $\alpha(3/2)$ as $\alpha(\text{FB})$

Organofluorine Selectivity indicator-2: $\alpha(6/4)$ as $\alpha(\text{TFT})$

➤ PFP-R showed high organofluorine selectivity [$\alpha(\text{FB})$, $\alpha(\text{TFT})$], comparable to PFP, with minimal batch variation.

With high affinity for organofluorine compounds, PFP-R is well suited for comprehensive **PFAS analysis**. Benefits include suppression of PFBA peak broadening, avoidance of matrix interferences by temperature screening, and improved separation of branched / linear PFAS.

Specifications

Stationary Phase	Particle size	Surface area	Pore diameter	Carbon content	End-capping	Available pH range	USP Code
SunBridge PFP-R	3 μm	190 m^2/g	15 nm	14%	Yes	2 - 8	L43

Ordering information

Packing	Inner diameter (mm)	2.1	2.1 [PS inert*]	3.0	4.6
	Length (mm)	Part number	Part number	Part number	Part number
SunBridge PFP-R 3 μm	50	JV2941	JV294PS	JV2341	JV2441
	100	JV2961	JV296PS	JV2361	JV2461
	150	JV2971	JV297PS	JV2371	JV2471
	250	JV2981	—	JV2381	JV2481

* PS inert is an inert column hardware option with a special surface that reduces metal adsorption (standard type: stainless steel).

Related products	Product name	Part number
Delay column for PFAS**	SunBridge Delay C18 3 μm 3 mm I.D. x 50 mm L	JB2341NR

** For PFAS LC/MS analysis, installed just after the pump to delay system-derived PFAS background peaks from the target analyte

Manufacturer

ChromaNik Technologies Inc.

6-3-1 Namiyoke, Minato-ku, Osaka, 552-0001 Japan

TEL: +81-6-6581-0885 FAX: +81-6-6581-0890

E-mail: info@chromanik.co.jp

URL: http://chromanik.co.jp